

Application:

Grouping ABO-RH1
Phenotyping RH-KEL1
Reverse typing

Internal Quality Control composed
of 4 known phenotyped samples



DESCRIPTION

DESIGNATION	REFERENCE	DESCRIPTION	STORAGE	STABILITY
CQI ABO-D - Rh.K	CQI02	4 tubes of whole blood phenotyped samples	2°C - 8°C	Until expiry date

ABO

1 sample
A1

1 sample
A2

1 sample
B

1 sample
O

RH

1 sample RH:
-1, -2, -3, 4, 5

1 sample RH:
1, 2, -3, -4, 5

1 sample RH:
1, -2, 3, 4, -5

1 sample RH:
1, 2, 3, 4, 5

KELL

At least
1 sample KEL:
1, 2 (K+ k+)

At least
1 sample KEL:
-1, 2 (K- k+)



Regulatory
CE marked



Lead time
6 weeks



Interpretation

